

User Manual for pager/responder

Welcome to choose the "Help" smart home series product's pager/ responder. The product consists of a pager, a responder, and a Wi-Fi LAN. Please read this manual carefully before using the product

Target audience for the product: elderly people or those in need of special care

Function

Users can send help requests via voice. The requests will be sent to the responder through the Wi-Fi local area network. After receiving the requests, the caregiver, that is, the person providing help, can respond via voice to inform the requester that the requests have been received.。

Installation guide

Please complete the system setup in the following order. **Important reminder: Please keep the product away from sources of fire, heat, water, and gas, and ensure that it will not cause harm to children!**

1. Install the mobile app. For specific instructions, please refer to the "Installation of mobile app" and "Mobile "app authorization" section later; Pay attention to the following key points: Authorize positioning function on your phone before installation, and authorize the app's location service permissions during installation; confirm the authorize app positioning function again after installation. The first time you run the app, a privacy policy window will appear. Please click "Agree" to enter the software.

2. Add the pager/responder master unit(There can only be one main unit)

Firstly, ensure that the phone is connected to the home's wireless network; please pay attention to the product label, select the master unit and use the adapter provided with it to power on. When you hear the voice prompt of "Distribution Network", click the three horizontal buttons in the upper left corner of the app page, and click "Add Pager" in the side menu. When the device "Carer" is connected, click "Connect". After a while, a prompt box "Add Completed" will pop up, indicating that the caller's main unit has been added. When you hear 'New slave unit access', it means that the communication between the master and slave units is smooth and can be used normally.

3. Add the pager/responder slave unit(Up to four slave units can be configured)

After completing the addition of the pager/responder master unit, please add the pager/responder slave unit. Please pay attention to the product label, select the slave unit, and use the adapter provided with it to power on. When you hear the voice prompt of "Distribution Network", click the three horizontal buttons in the upper left corner of the app page, and click "Add Pager" in the side menu. When the device "Carer" is connected, click "Connect" and wait for a prompt box of "Add Completed" to pop up, indicating that the operation is completed. When you hear 'Connect main unit', it means that the communication between this unit and the main unit is smooth and can be used normally

Usage

When in use, the pager should be placed in the room of the person in need, and the responder should be placed in the room of the caregiver (who provides assistance). Both the master and slave units can be used as pager or responder.

Send a request for help when the user needs help, simply call out clearly 'Help', and when you hear the response 'I'm here, master', please say 'come and help me'. Under normal circumstances, Seeking help information is sent in real-time to the responder. Similar distress commands include common needs such as "I'm thirsty", "I'm hungry", "I need to go to the bathroom", etc.

Response to Help Request When the responder receives a request for help, the caregiver can answer "I know" or "received" or "here" to notify the help seeker that the request has been confirmed.

Volume adjusting you can first call 'someone here', and after receiving a response from 'here, master', say 'increase volume' or 'decrease volume' to complete the volume adjustment. Even if the power is turned off, the adjusted volume will remain unchanged until the next volume adjustment.

Installation of mobile app

This section and the following mobile authorization operations may have different expressions depending on the brand and model, and the corresponding operations should be completed according to the specific situation.

Before installing the mobile app, please Authorized positioning function on your phone. Press and hold the top of the phone screen or swipe down from the top right corner to select location services and turn them on.

Enter the mobile app store, search for 'Carer', or use your phone's browser to scan the QR code in the manual and install the app. During the installation process, you will receive requests for location services and message notification permissions. Please select 'Allow'. Complete the installation according to the installation prompts.

Mobile app authorization

Authorize positioning function

Follow the steps below on your phone: Settings - Apps - App Management - Select "Carer" - Permissions - Location - Allow When Using

Jilin Advantech Computer Technology Co., Ltd